

# "I've stopped falling!"

## Finally! a walking aid for Parkinson's

"For me it is wonderful. It keeps me from falling and enables me to walk, and thus exercise my legs. I no longer have to have help to walk."  
— Mrs. L.D. McDade

### Designed for Parkinson's Disease!

#### • Stability

We created the *U-Step Walking Stabilizer* to increase your independence. The secret is in our patented U-shaped base. Its ultra-stable foundation braces you in every direction. It is not like pushing a walker. Instead, the *U-Step* surrounds you and moves with you. You will feel as stable as you would feel holding onto another person's arm.

#### • Safety

The innovative braking system is easy to use and puts you in complete control. The *U-Step* will not roll unless you are ready to walk. When you lightly squeeze a hand brake, the unit will roll with you. Once you release the hand brake the unit will stop immediately. This feature is particularly helpful when standing up from a chair because, the unit will not roll away from you.

- Rolls easily
- Improves posture
- Improves walking stride
- Rolling resistance control
- Folds for transport
- Easy turning

"I recommend the *U-Step*. It is useful for patients with balance problems, because of its wide base."

— Abraham Lieberman, M.D., Medical Director  
National Parkinson Foundation, Miami, FL



Optional feature:  
red button for operating  
LaserLight™

### Laser Light™ helps prevent freezing

If you suffer from freezing episodes, *U-Step* has a solution. Introducing *LaserLight™* an exclusive optional feature of the *U-Step*. Simply press the red button on the handlebar and a bright red laser line is projected on the floor for you to step over. See *projected line illustrated on back page*. It is amazing! Our laser offers an entirely safe, obstacle-free, visual cue that helps you break the "freezing episode," and walk normally with increased stride length.

"My father shows off the laser to all his cohorts with Parkinson's. He says that if you can throw down a pencil and step over it, the laser will work for you. He does not use the urinal at the side of his bed in the middle of the night any more. With the *U-Step* and the laser he can make it to the washroom."

— Dr. Carlo Klott, M.D., Austin, TX

Note: None of the physicians or therapists quoted in this material have any financial interest in selling the *U-Step*.



**Stop falling! Order your *U-Step* or *LaserCane* now!**

***LaserCane* is ideal for those seeking minor walking support**

