

AQUATEC OCEAN THE NOVELTY

OCEAN basic model



- ... Stable, stainless steel frame
- ... Individually adjustable seat height (simple to adjust, no tools required)
- ... Washable backrest with quick-release buckles – can be adjusted to fit the user's needs
- ... Comfortable due to the ergonomically designed seat
- ... Removable seat – easy to clean
- ... Height-adjustable, swivel-mounted footrests for easy transfer and adjustment
- ... Structured fold-away footrests with adjustable heel loops provide a safe hold
- ... Four braked castors for reliable stopping
- ... Easily manoeuvrable due to smooth castors
- ... Ergonomic fold-away armrests
- ... Knock-down concept allows easy transport and use on travels
- ... Space-saving packing size
- ... Comprehensive accessories for specific requirements
- ... Smooth surfaces and flowing contours for easy cleaning



AQUATEC OCEAN **AQUATEC OCEAN**
24"-WHEELS

Order number	7.10.900	7.12.900
Dimensions		
Height to armrests (A)	680-810 mm	680-810 mm
Seat height (B)	477-600 mm	477-600 mm
Seat depth (C)	450 mm	450 mm
Overall width (D)	560 mm	675 mm
Overall length (E)	934 mm	1060 mm
Length without footrests (F)	730 mm	856 mm
Footrest depth (G)	154 mm	154 mm
Overall height (H)	964-1094 mm	964-1094 mm
Height seat to lockable footrest (I)	280-530 mm	280-530 mm
Height ground to lower seat edge:		
With holder for sanitary pan (J)	390-520 mm	390-520 mm
Without holder for sanitary pan (J)	410-540 mm	410-540 mm
Width between armrests (K)	459 mm	459 mm
Height adjustment steps	25 mm	25 mm
Wheel diameter	127 mm	127 mm
Wheel diameter of 24" self-propelling wheels	–	610 mm
Others		
Weight	14.5 kg	19.5 kg
Load capacity	130 kg	130 kg
Colour	white/blue	white/blue

Heel loops are included in delivery of the **AQUATEC OCEAN**. Sanitary pan holder and sanitary pan are not included. Other accessories can be found on pages 15 and 16.



OCEAN with 24" self-propelling wheels

