

Premium Features are Standard on Breezy Ultra⁴




Quick-Release Wheels
Make loading in/out of vehicles easy.



Angle-Adjustable Front Casters
Accommodate angled seating surfaces without moving the caster housings out of square.



Adjustable Axle Positions
Three seat height positions and the ability to add a fixed tilt allow for greater flexibility.



Wraparound Armrest
Provides a narrower overall chair width.



Quick-Adjust Footrests
Adjust height without tools using spring-loaded detent buttons. Heel loops adjust with Velcro[®]



Height-Adjustable Back with Pocket
Store personal items/documents in pocket with closing flap. Adjust height easily without tools.

Key Options




Angle-Adjustable Back
Back adjusts from 6° anterior to 24° posterior in 6° increments for better positioning. Height adjusts to 17", 18" or 19".



Height-Adjustable Arms
4" vertical adjustment allows for use of cushions.

Colors



Actual colors may appear different.

Specifications

Frame Weight*	30.9 lbs (14.02 kg)	
Frame Material	T6 6061 Aluminum	
Seat Width & Depth	16"x16"	16"x18"
	18"x16"	18"x18"
	20"x16"	20"x18"
Back Height	Adjustable	17", 18", 19"
Seat Height	15.5", 16.5", 17.5", 18.5", 19.5"	
Upholstery	Padded Seat	Non-Padded Seat
	Nylon Padded Back w/wo Pocket	
Front Casters	6" Polyurethane	
	6" Solid	
	8" Polyurethane	
	8" Solid	
Rear Wheels	20" Full Profile Polyurethane	
	with Aluminum Handrims 24" Low Profile Polyurethane	
Axles	Quick Release	
Armrest	Flip-Back, Desk or Full	
	Flip-Back Height Adjustable, Desk or Full	
Front Hangers	Swing-Away with Heel Loops	
	Elevating Legrest	
Maximum User Weight	250 lbs (113.4 kg)	
Color Options	Candy Blue, Black, Forest Green, Silver-Texture, Silver Polish	
Accessories	Seat Belt, Auto Style	BZU9001
	IV Pole	BZU9003
	Tray Table, Full Arm	BZU9002
	Tray Table, Desk Arm	BZU9009
	Oxygen Tank Holder	BZU9004
	Crutch Holder	BZU9005
HCPCs Code	K0004 (meets requirements for)	

* 18"x16" Configuration without front hangers