

# Premium Hi-Lo Adjustable Bed with Trendelenburg option

**We know the value of a  
Good Night's Sleep**  
*and the importance of caring  
for loved ones*

**15**  
Year  
Warranty  
On Steel Frame

Liberty HealthCare's PREMIUM HI-LO with TRENDELBURG option Adjustable Bed is a highly specialised bed designed for high care patients with special needs. It has a higher raising capacity than the HOMECARE HI-LO and the Trendelenburg function tilts the entire bed forward or backwards to assist with patient care and treatment. It features Linak motors and a high quality steel frame to give you many years of trouble free use.

Featuring:

- 550mm lift capability of Hi-Lo Function (vertical raising and lowering of the bed)
- Trendelenburg option – Forward and Reverse tilt
- Direct Drive and Twin Drive Motors
- Soft Touch Corded Handset



Liberty HealthCare PREMIUM HI-LO with TRENDELBURG beds can be custom made to suit the user and carer's specifications. A large range of accessories are available to make this very functional bed look right at home in your bedroom, including bed rails and head & foot boards. Ask about our companion sets available throughout our range.

Small Single 2'6"	Single 3'0"	King Single 3'6"	Double 4'6"	Queen 5'0"
-------------------------	----------------	------------------------	----------------	---------------

Split Queen 2 x 2'6"	Split King 2 x 3'0"
----------------------------	---------------------------

**Better Sleep = Better Health = Better Life**

