

Shoprider GK Series

ASK FOR THE BRAND SHOPRIDER



GK10 CROSSOVER

GK9 LITTLE RIPPER

SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®

www.shoprider.com.au

mobility

GK Series

The SHOPRIDER™ **GK Series** of scooters are designed to be fully portable and offer exceptional functionality to the customer. Each of the **GK Series** scooters offer many of the same features from the larger models all with the advantage of being portable. Every model in the **GK Series** has the same easy pull apart functionality, fully rotating lockable seat, adjustable steering tiller and fold up arm rests.

The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	GK9	GK10	FEATURES	GK9	GK10
Length	1050mm	1285mm	Mid Back Captain's Seat	N	Y
Width	560mm	550mm	Mid Back Seat	Y	N
Height	880mm	960mm	Adjustable Steering Tiller	Y	Y
*Safe Climbing Angle	9 Degrees	9 Degrees	Adjustable Width Arm Rests	Y	Y
*Max Climbing Angle	12 Degrees	9 Degrees	Fold up Arm Rests	Y	Y
*Maximum Load Single user	136Kgs	136Kgs	Sliding Seat	N	Y
Charger (Del with min)	4Amp	4Amp	Rotating Seat	Y	Y
Motor Size	0.75Hp	0.8Hp	Front Shopping Basket	Y	Y
Maximum Speed	6Kph	8Kph	Variable seat height	Y	Y
Battery capacity min	2 x 21Amp	2 x 36 Amp	Horn	Y	Y
Maximum range #	up to 20km	up to 28km	Battery Level Warning Light	N	Y
Wheel Size	220mm	260mm	Battery Level Gauge	Y	N
Heaviest Part	16.5kg	21kg	Front Lights	N	Y
Weight (w/o Batteries)	41.5kg	49kg	Split Battery Box	Y	Y
			Pneumatic Tyres	N	Y
			Solid Tyres	Y	N
			Anti Tips Wheels	Y	Y
			Comprehensive Warranty	Y	Y


Shoprider™ is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

WARNING: Do not exceed 10kph on public roads or Paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



*(Single user)

GK9	GK10
12° (20%) 70kg	
9° (15%) 90kg	9° (15%) 60kg
6° (10%) 113kg	6° (10%) 85kg
0° (0%) 136kg	0° (0%) 113kg

SHOPRIDER®

www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au