



FEATURES OF THE BARIATRIC GOWNS

- Gowns are Non see-through, comfortable and breathable with long sleeves.
- The plastic buttons on the sleeves open up both sides, enabling easy access for procedures but also providing privacy for the patient.
- Designed using a three layer dark blue SMS fabric, 50gsm which is lightweight.
- Includes a secure tie, is ideal for Hospital Patients, Day Procedures, GP clinics, Radiology or any setting where a patient is required to wear a disposable gown.
- Versatile sizing to accommodate all body shapes with a waistline up to 2.5m.

SIZES AVAILABLE

When ordering please quote the following
 (Minimum quantity 1 box):

DESCRIPTION	SIZE/COLOUR	QTY/BOX
Disposable Bariatric Gown: XXXL	L 120cm x W 200cm/Dk Blue	50
Disposable Bariatric Gown: XXXXL	L 120cm x W 240cm/Dk Blue	50

GP Medical

Unit 2/2 Enterprise Cres
 Malaga WA 6090

Tel: (08) 9249 9333
www.gp-medical.com.au