

# IndeeLift®

The Human Floor Lift

For Professional Healthcare

*Acute Hospitals  
Skilled Nursing  
Assisted Living*

The People Picker Upper!

HFL-550P



*"IndeeLift is the tool you need to safely get your residents off the ground after a fall."*





# The People Picker Upper!

## HFL-550P

### Purpose

The HFL-550P models is a portable human floor lifts designed for use in professional healthcare environments to safely lift patients and residents who have fallen and are unable to get up from the floor. These sturdy, stable lifts utilize medically certified components with a working load of 250kg

### Operation

After a fall, the healthcare professional wheels the HFL to the fallen person and assists them onto the seat. With the simple push of a button on the wired remote, the fallen person is safely and easily raised from the floor.

### Summary of Features

- Lifts a person weighing up to 250kg to a height of 53.4cm above the floor.
- No complex cradles or harness to configure.
- Easy maneuverability for flexible positioning in tight spaces including ADA-compliant bathrooms.
- Large ramp-mount seat for easy slide-on access.
- Swinging/removable Rise Assist Handles for easy transfer to wheelchair.
- Powder coated for easy cleaning and decontamination.
- Powered by a battery pack with built-in charger.
- Audible indicator when charging is required.
- Charges from any standard AC power outlet.



### Product Specifications

#### Physical Specifications

Safe Working Load.....  
 Maximum Overall Length.....  
 Maximum Overall Height.....  
 Minimum Overall Height.....  
 Maximum Seat Height.....  
 Turning Radius.....  
 External Width .....  
 Overall Height of Legs.....  
 Rear Casters.....  
 Power Pack.....  
 Total Weight.....

#### HFL-550

550lbs .....250kg  
 29"..... 735mm  
 64".... 1626mm  
 45".... 1143mm  
 21".....534mm  
 32".....813mm  
 26.25"...667mm  
 2.5".....64mm  
 5".....127mm  
 6.6lbs.....3kg  
 108 lbs.....49kg

All measurements are within a +5/-5% tolerance.

#### Electric Shock Protection

Charger – Class II  
 Lift – Internal power source

#### Degree of Shock Protection

Charger – Type B  
 Lift – Type B  
 Operating environment: >+5°C <+40°C

#### Electrical Specifications

Batteries – 2x12 volt Rechargeable sealed lead acid type  
 Battery Capacity – 3.2 Amp Hours  
 Charger Rated Input – 100/240VAC 50/60Hz  
 Charger Rated Output – 27.4/29 VDC@0.8A

**Tel: 08 9249 9333**  
**Email: info@gp-medical.com.au**

