

HEIGHT ADJUSTABLE
180mm



**NEW
PRODUCT
RELEASE!**

**REMOVABLE
ARMRESTS / WINGS**
For SIDE TRANSFERS

**REMOVABLE
LITHIUM-ION
BATTERY**
No Cables

MOBILE
Ø75mm Castors



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ULTRACARE MOBILE LIFT CHAIR

LC0901

FEATURES

- Lift / Rest / Recline armchair with height adjustment.
- Mobile lift chair with removable Lithium-ion battery.
- Ø75mm castors with push handle at rear.
- Handset controls with large buttons.
- Removable armrests and wings for side transfers.
- Incontinent PU backrest / seat / legrest coverings
- Vinyl side coverings with storage pockets either side.
- Added comfort of a pocket spring seat base.
- Quality German actuators.
- Warranty: Frame: 5 years / Battery: 1 year / Electrics: 2 years / Fabric: 2 years

SPECIFICATIONS All dimensions shown are nominal

Overall: **760mm (W) x 850mm (D) x 1140mm (H)**

Seat: **540mm (W) x 510mm (D) x 530 - 710mm (H)**

Armrest: **100mm (W) x 450mm (L) x 180mm (H)**

Backrest: **540mm int. (W) x 700mm (H): 110 - 145°** (Relevant to seat)

Legrest: **540mm (W) x 450mm (L): 90 - 15°** (Relevant to seat)

Castors: **Ø75mm - Swivel front / Braked rear**

Distance from wall: **600mm**

Product Weight: **65kg**

Maximum User Weight: **160kg**

Power Source: **Removable rechargeable Lithium-ion battery pack**
25.2V - 2550mA - 65Wh



G&P MEDICAL ENTERPRISES

P: 08 9249 9333

E: sales@gp-medical.com.au